

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

MBD2 (Methyl-CpG-Binding Protein 2) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01302

---

第三版

**描述:** Methylation at CpG dinucleotides in genomic DNA is a fundamental epigenetic mechanism of gene expression control in vertebrates. Proteins with a methyl-CpG-binding domain (MBD) can bind to single methylated CpGs and most of them are involved in transcription control.

**抗原:** Synthetic peptide of part of C terminal of human MBD2 (TEEMDIEMDSGDEA)

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Not cross-reacts with MBD3 (by western blotting)