

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MDA5 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01304

第三版

描述:A full-length cDNA encoding MDA5 was isolated by screening a melanoma differentiation-subtracted library and a placenta cDNA library, also known as Ifih1 (Interferon induced with helicase C domain protein 1).

抗原:Synthetic peptide of part of N terminal of human MDA5 (TIEDRNRIAAAENNGN)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Reacts with human and mouse MDA5