



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Rab27A/B Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01305

第三版

描述: Small GTPase Rab is a large family of membrane trafficking proteins that are conserved in all eukaryotic cells. More than 60 Rab isoforms have been reported in mice and humans, and they are believed to regulate various steps (or various types) of organelle transport.

抗原: Recombinant Rab27B protein

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Cross-react with mouse and rat Rab27A/B. Not cross-react with the Rab3 subfamily