
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MIP-2 (C) Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01307

第三版

描述: Cytokine-induced neutrophil chemo attractant 1 (CINC-1) was originally purified from media conditioned by IL-1 stimulated rat kidney epithelioid cells (NRK-52E). Amino acid sequence that encodes for rat CINC-1 was identified in 1989 by Watanabe ' s group at Toyama Medical and Pharmaceutical University.

抗原: Synthetic peptide for C-Terminal of Mouse MIP-2

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Mouse MIP-2 specific, Non-cross react with Mouse KC.