

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

mLST8 (L17) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01310

第三版

**描述:**mLST8 (mammalian lethal with sec-13), an about 39 kDa protein which has WD40 domains, has been identified as a configuration factor of mTORC1/mTOR complex 1 and then found to be also a configuration factor of mTORC2.

**抗原:**Synthetic peptide of N-terminal of Human mLST8 (TAGYDHTVRFWQAHSGLIC)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**