



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse IgG (H+L) Goat IgG Affinity Purify

多克隆抗体

产品货号: PA01312

第三版

描述:Anti-Mouse IgG (H+L) Goat IgG Affinity Purify, 1 mg. Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃, purified with antigen.

抗原:Mouse IgG (H+L)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:The antibody is absorbed by rat IgG, goat IgG, rabbit IgG and human IgG. Cross-reactivity to goat IgG, rabbit IgG and human IgG is less than 0.1 %. Cross-reactivity to rat IgG is about 0.5 %.