

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Raptor (R1) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01314

第三版

**描述:**Raptor (regulatory associated protein of mTOR) is an about 150 kDa protein which has HEAT repeats and WD40 repeats, and it is identified as an mTOR binding protein. It is conserved over species and widely expressed in a variety of cells.

**抗原:**Synthetic peptide of the N terminal part of Human Raptor (MESEMLQSPLLGLGEEDEA)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**