

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MRS (Methionyl-tRNA synthetase )

多克隆抗体

产品货号: PA01315

第三版

**描述:**MRS is a component of multi-ARS complex and catalyzes ligation of methionine with methionyl-tRNA. MRS is involved in the rRNA synthesis in nucleoli.

**抗原:**Denatured full length human recombinant MRS.

**配方:**Anti-methionyl-tRNA synthetase rabbit serum containing 0.05% sodium azide for preservative.

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**This antibody reacts specifically with 100kDa methionyl tRNA synthetase protein.