
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

N-Cadherin (760) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01321

第三版

描述: Cadherins are a family of Ca²⁺-dependent intercellular adhesion molecules. They play an important role in cell-cell interaction, histogenesis and cellular transformation, and association with the actin cytoskeleton regulates their function by a complex of cytoplasmic proteins called the catenins (, ,).

抗原: Synthetic peptide of aa 746-761 (KRRDKERQAKQLLI) of human N-Cadherin

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: