
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-RET Finger Protein (RFP) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01322

第三版

描述:RET finger protein (RFP) belongs to the large B-box RING finger protein family and is known to become oncogenic by fusion with RET tyrosine kinase. Although RFP is reported to be a nuclear protein that is present in the nuclear matrix, its function is largely unknown.

抗原:Synthetic peptide for C-terminal portion (IDGFSGHVGNHGHSMETSP) of RFP

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:The antigen peptides is selected a common region between human and mouse