
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

N-Cadherin (YS) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01323

第三版

描述: Cadherins are a family of Ca²⁺-dependent intercellular adhesion molecules. They play an important role in cell-cell interaction, histogenesis and cellular transformation, and association with the actin cytoskeleton regulates their function by a complex of cytoplasmic proteins called the catenins (, ,).

抗原: Recombinant human N-Cadherin extracellular subdomain 1-3

配方: Lyophilized product from PBS solvent (BSA and NaN₃ are not containing)

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Non-cross reacts with another cadherin ' s transfectants