
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Rheb (R165) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01328

第三版

描述:Rheb (ras homolog enriched in brain) is a 20 kDa protein identified as a Ras homolog expressed in the rat hippocampus, but it widely exists in a variety of cells. It is also reported that GTP-binding form is predominant in Rheb, unlike other low molecular weight GTP-binding proteins.

抗原:Synthetic peptide of the C-terminal part of Human Rheb (LEAEKMDGAASQGKSS)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: