
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Rictor (R1291) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01330

第三版

描述:Rictor (rapamycin-insensitive companion of mTOR/ mAVO3) is about 200 kDa and identified as an mTOR-binding protein. Rictor is a factor that constitutes a complex with mLST8 and Sin 1 which is called mTORC2 (mTOR complex 2), and it is contrasted with raptor in mTORC1.

抗原:Synthetic peptide of a part of Human Rictor (GSSHTLPRRAQSLKA)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: