
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

RON (G1385) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01334

第三版

描述:The c-MET gene encodes the tyrosine kinase p190MET, the receptor for a molecule known as Scatter Factor (SF) or Hepatocyte Growth Factor (HGF). This molecule has biological activities on epithelial sheets, including mitogenesis, cell-cell dissociation, stimulation of migration into the extracellular matrix, induction of cell polarization and branched tubulogenesis.

抗原:Synthetic peptide for a part of C terminal of human RON

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Confirmed its specific reaction with RON/BaF3 transfectant