

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Nephrin (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01335

第三版

描述: Nephrin has been identified as a disease causing molecule responsible for congenital nephrotic syndrome of Finnish type, in which the proteinuria starts in utero or immediately after birth (ref. 1). It is an immunoglobulin-like transmembrane adhesion protein having eight Ig motifs. Nephrin forms a structural backbone of the slit diaphragm through a head-to-head, trans-homo-interaction at the N-terminal ectodomains.

抗原: Synthetic peptide of the C-terminal part of Human Nephrin

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Cross-reacts with rat.