



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Nephrin (N) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01337

第三版

描述: Nephrin has been identified as a disease causing molecule responsible for congenital nephrotic syndrome of Finnish type, in which the proteinuria starts in utero or immediately after birth (ref. 1).

抗原: Recombinant protein of Human Nephrin

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Please see datasheet.