

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RRS (Arginyl-tRNA synthetase)

多克隆抗体

产品货号: PA01338

第三版

描述:RRS is a cofactor of multi-ARS complex and mediates ligation of arginine with tRNA^{Arg} using ATP and is associated with AIMP1 in ARS complex.

抗原:Native 1-72aa recombinant human RRS.

配方:Anti-arginyl-tRNA synthetase rabbit serum containing 0.05% sodium azide for preservative.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:This antibody reacts specifically with 74kDa arginyl tRNA synthetase and immunoprecipitates with AIMP1 (p43).