



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Nestin (N1602) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01340

第三版

**描述:** Nestin is one of the intermediate filaments abundantly produced in the developing central nervous system and sometimes in the embryonic stage. Nestin is also reportedly detected in gliomas/glioblastomas. Thus, the antibody against nestin can be widely used as a marker for proliferating neural progenitor cells in the nervous system.

**抗原:** Synthetic peptide for C Terminus (KFTQREGDRESWSSGED) of human Nestin

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Confirmed by western blotting.