



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-S-100 Protein Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01342

第三版

描述:S100 proteins (16 members) show a very divergent pattern of cell- and tissue-specific expression, of subcellular localizations and relocations, of post-translational modifications, and of affinities for Ca^{2+} , Zn^{2+} , and Cu^{2+} , consistent with their pleiotropic intra- and extracellular functions.

抗原:Bovine Brain S-100 Protein

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN_3

如何使用:加1ml超纯水重溶

稳定性: -20°C 保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 $\mu\text{g/ml}$.

特异性:React to α and β chain of Bovine and Human S-100 protein