
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-SCF (K165) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01347

第三版

描述:Stem cell factor (SCF) and its receptor KIT play an important role in various biologic phases, such as hematopoiesis, reproduction, and regeneration. Structurally, c-Kit contains five immunoglobulin-like domains extracellularly and a catalytic domain divided into two regions by a 77 amino acid insert intracellularly.

抗原:Synthetic peptide for a part of Human SCF

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Crossreact to Mouse and Rat SCF. Positive to recombinant SCF(E.Coli). Not tested about reaction to native membrane bound type