

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-SCF (K236) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01349

第三版

**描述:**Stem cell factor (SCF) and its receptor KIT play an important role in various biologic phases, such as hematopoiesis, reproduction, and regeneration. Structurally, c-Kit contains five immunoglobulin-like domains extracellularly and a catalytic domain divided into two regions by a 77 amino acid insert intracellularly.

**抗原:**Synthetic peptide for a part of Human SCF

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Crossreact to Mouse SCF. Positive to recombinant SCF