

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Niban (53) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01350

---

第三版

**描述:** Niban has been cloned from cell lines derived from Eker rat by AFLP-Differential Display method as a novel gene which is increased in the early stage of renal carcinogenesis (ref. 6). Then, it has been reported that Niban protein is expressed in various cancer tissues including thyroid cancer. In research of recent years, it has been reported that Niban is an ER stress-inducible protein and it regulates phosphorylation of eIF2 and S6K1/4E-BP1 (ref. 3).

**抗原:** Synthetic peptide of a part of rat Niban (TDLLPKLKGKKNDRKRA).

**配方:** Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Please see datasheet.