

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-SEP1 Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01352

第三版

描述:SEP1 has been reported by group of Takahashi, as a gene enhancing its expression by stimulation by a neurotrophic factor, GDNF. Though SEP1 is found in every organ, it is located at neurite of neuroblastoma cells of which differentiation is induced by GDNF, and SEP1 expression has been confirmed in many of neurons in an organism. That indicates a possibility that SEP1 relates to neural functions.

抗原:Synthetic peptide of the part of Human SEP1

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:React with human and mouse SEP1