
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Septin 4 (N) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01359

第三版

描述: Among the members of the mammalian septin family, septin 4 (H5) is unique in terms of its primary structure and expression pattern. While its expression has also been demonstrated in the heart, liver and testis, it is most prominently expressed in some neurons and glial cells.

抗原: The N terminal part of synthetic peptide of Human Septin 4 (GKLDPYDSSEDDKEYVG)

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Cross reacts with mouse and rat.