
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Septin 7 (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01361

第三版

描述:A cytoskeletal or scaffold protein family named “ septin ” has been shown to be involved in diverse phenomena in the living body, including cell division. The term “ septin ” is derived from “ septation, ” which means division.

抗原:Synthetic peptide of the C terminal part of Human Septin 7 (NSSRTLEKNKKKGKIF)

配方:Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Cross reacts with mouse and rat.