

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Single Stranded DNA (ssDNA) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01363

---

第三版

**描述:** Apoptosis, or programmed cell death, is a highly conserved biochemical mechanism that allows cells to die in a controlled and organized manner. This death process is essential for normal cellular differentiation and tissue homeostasis within multi-cellular organisms.

**抗原:** Bovine fragment single stranded DNA

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Can reacts with random sequenced single DNA of all species not only of human.