
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Sir2/SIRT1 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01364

第三版

描述: Imai et al. have found that *S. cerevisiae* Sir2p and its higher eukaryotic orthologs are novel nicotinamide adenine dinucleotide (NAD)-dependent histone/protein deacetylases (1). Sir2 proteins promote longevity in yeast and in *C. elegans* by silencing rDNA and regulating the insulin/IGF-1 signaling pathway, respectively. The

抗原: Synthetic peptides for part of C Terminal of human Sir2/SIRT1 (LEDEPDVPERAGG)

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Confirmed by western blotting.