

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Smad3C (Ser 423/425 Phosphorylated) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01367

第三版

**描述:**Smads are critical mediators conveying signals from the TGF- $\beta$  superfamily members to the nucleus (ref. 1-3). Activated TGF- $\beta$  type I receptor (T $\beta$ RI) and c-Jun N-terminal kinase (JNK) differentially phosphorylate Smad3 to convert two distinctive phosphoisoforms.

**抗原:**Synthetic peptide of phosphorylated Smad3C (Ser 423/425)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**