
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Olig2 Anti-Human Rabbit IgG Affinity purify

多克隆抗体

产品货号: PA01371

第三版

描述:Oligodendrocytes are defined as the cells that exist in tissue stroma of central nervous system and play a important role in maintaining function of neuron cell by making myelin wrapping around a nerve cell axon.

抗原:Synthetic peptide in portion of C-terminus of Human Olig2

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Cross-reacts with Mouse and Rat