



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Osteopontin (O-17) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01375

第三版

**描述:**Osteopontin (OPN) is a secreted glycoprotein that was originally isolated from bone. Its molecular weights have been reported in the range of 66 kDa to 44 kDa depending on glycosylation and phosphorylation.

**抗原:**Synthetic peptide for N-terminal of Human Osteopontin

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**React to both recombinant and native human Osteopontin. Can be react with thrombin digested form of human Osteopontin. Non-Cross react to mouse Osteopontin