
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Osteopontin (O-17) Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01377

第三版

描述:Osteopontin (OPN) is a secreted glycoprotein that was originally isolated from bone. Its molecular weights have been reported in the range of 66 kDa to 44 kDa depending on glycosylation and phosphorylation.

抗原:Synthetic peptide for N-terminal of Mouse Osteopontin

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:React to both recombinant and native mouse Osteopontin. Cross react to Rat Osteopontin and non-cross react to Human Osteopontin