
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Osteopontin (O-17) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01382

第三版

描述: Osteopontin (OPN) is a secreted glycoprotein that was originally isolated from bone. Its molecular weights have been reported in the range of 66 kDa to 44 kDa depending on glycosylation and phosphorylation.

抗原: Synthetic peptide for N-terminal of Rat Osteopontin

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: React to both recombinant and native rat Osteopontin. Detectable to Rat Osteopontin Thrombin digested form. Non-cross react to Human and mouse Osteopontin