

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Podocin Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01398

---

第三版

**描述:**Podocin has been identified as a disease causing gene responsible for familial steroid-resistant nephrotic syndrome (focal segmental glomerulosclerosis) (ref. 1). It is an integral membrane protein having a hairpin-loop-like intramembrane topology and is exclusively expressed at the base of the podocyte slit diaphragm (ref. 2).

**抗原:**Synthetic peptide of a part of Human Podocin

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Cross-reacts with rat and mouse.