

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Podoplanin Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01401

第三版

描述:Podoplanin belongs to a type-I transmembrane protein which is a sialoglycoprotein with some O-glycosylation sites. It is known as a lymphatic endothelium-specific marker and believed to have a platelet-aggregating activity and to promote hematogenous metastasis of cancer cells. Podoplanin is thought to be expressed on the shape alteration of squamous cells and to relate to transition of cytoskeletal actin structure. Also, its activity which decreases the adhesion ability of cells and promotes their metastasis has drawn researchers' attention. Additionally, podoplanin is reported to be expressed in human tumors such as brain tumor, oral cavity cancer, head and neck cancer, esophageal cancer, lung cancer, testicular tumor, cervical cancer and malignant mesothelioma.

抗原:Synthetic peptide of the part of Human Podoplanin

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: