



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PRAS40 (P238) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01404

第三版

描述:PRAS40 (proline-rich Akt substrate of 40 kDa) is a protein with a prolin-rich domain and appears at about 40 kDa by SDS-PAGE. Initially it was identified as a substrate of PKB/Akt, it was later shown that it is also phosphorylated by mTORC1 (a complex including scaffold protein which is necessary for mTOR and phosphorylation of substrate by mTOR).

抗原:Synthetic peptide of the C terminal part of Human PRAS40 (DLPRPRLNTSDFQKLKRY)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Please see datasheet.