

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Acetyl-HSPA1A/B (Lys246) Polyclonal Antibody

多克隆抗体

产品货号: PA01446

第三版

**描述:**In cooperation with other chaperones, Hsp70s stabilize preexistent proteins against aggregation and mediate the folding of newly translated polypeptides in the cytosol as well as within organelles. These chaperones participate in all these processes through their ability to recognize nonnative conformations of other proteins. They bind extended peptide segments with a net hydrophobic character exposed by polypeptides during translation and membrane translocation, or following stress-induced damage. In case of rotavirus A infection, serves as a post-attachment receptor for the virus to facilitate entry into the cell.

**抗原:**Synthesized peptide derived from the Internal region of human HSP70 around the acetylation site of K246.

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**