

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat/Monkey Phospho-CDK1/2/3 (Thr14) Polyclonal Antibody

多克隆抗体

产品货号: PA01524

第三版

描述: Plays a key role in the control of the eukaryotic cell cycle. It is required in higher cells for entry into S-phase and mitosis. p34 is a component of the kinase complex that phosphorylates the repetitive C-terminus of RNA polymerase II.

抗原: Synthesized peptide derived from human Cdk1/2/3 around the phosphorylation site of Thr14

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: breast cancer tissues, cancer cell lines and glioblastoma tissue.