
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Phospho-Cyclin E1 (Thr395) Polyclonal Antibody

多克隆抗体

产品货号: PA01539

第三版

描述:Essential for the control of the cell cycle at the G1/S (start) transition.

抗原:Synthesized peptide derived from human Cyclin E1 around the phosphorylation site of Thr395

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in testis and placenta. Low levels in bronchial epithelial cells.