

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse Phospho-CK-8 (Ser73) Polyclonal Antibody

多克隆抗体

产品货号: PA01540

---

第三版

**描述:** Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

**抗原:** Synthesized peptide derived from human Cytokeratin 8 around the phosphorylation site of Ser73

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Observed in muscle fibers accumulating in the costameres of myoplasm at the sarcolemma membrane in structures that contain dystrophin and spectrin. Expressed in gingival mucosa and hard palate of the oral cavity.