

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-MAP2K4 (Ser80) Polyclonal Antibody

多克隆抗体

产品货号: PA01593

第三版

描述:Dual specificity kinase that activates the JUN kinases MAPK8 (JNK1) and MAPK9 (JNK2) as well as MAPK14 (p38) but not MAPK1 (ERK2) or MAPK3 (ERK1).

抗原:Synthesized peptide derived from human MEK-4 around the phosphorylation site of Ser80

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Abundant expression is seen in the skeletal muscle. It is also widely expressed in other tissues.