

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat Phospho-MAP2K4 (Thr261) Polyclonal Antibody

多克隆抗体

产品货号: PA01594

---

第三版

**描述:** Dual specificity kinase that activates the JUN kinases MAPK8 (JNK1) and MAPK9 (JNK2) as well as MAPK14 (p38) but not MAPK1 (ERK2) or MAPK3 (ERK1).

**抗原:** Synthesized peptide derived from human MEK-4 around the phosphorylation site of Thr261

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Abundant expression is seen in the skeletal muscle. It is also widely expressed in other tissues.