

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Phospho-PDLIM5 (Tyr251) Polyclonal Antibody

多克隆抗体

产品货号: PA01620

第三版

描述: May play an important role in the heart development by scaffolding PKC to the Z-disk region. May play a role in the regulation of cardiomyocyte expansion. Overexpression promotes the development of heart hypertrophy. Contributes to the regulation of dendritic spine morphogenesis in neurons. May restrain postsynaptic growth of excitatory synapses.

抗原: Synthesized peptide derived from human PDLIM5 around the phosphorylation site of Tyr251

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Heart and skeletal muscle specific. Expression is commonly increased in the brain of patients with bipolar disorder, schizophrenia, and major depression.