

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-PRAG1 (Tyr413) Polyclonal Antibody

多克隆抗体

产品货号: PA01623

第三版

描述: This gene encodes an enzyme that belongs to the tyrosine protein kinase family. A similar protein in rat binds to Rho family GTPase 2 (Rnd2) and regulates neurite outgrowth via activation of Ras homolog gene family, member A (RhoA). PRAG1 (PEAK1 Related Kinase Activating Pseudokinase 1) is a Protein Coding gene. An important paralog of this gene is PEAK1.

抗原: Synthesized peptide derived from human Pragmin around the phosphorylation site of Y413.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: