

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat/Monkey Phospho-RAF1 (Ser259) Polyclonal Antibody

多克隆抗体

产品货号: PA01625

---

第三版

**描述:**Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. Part of the Ras-dependent signaling pathway from receptors to the nucleus. Protects cells from apoptosis mediated by STK3.

**抗原:**Synthesized peptide derived from human Raf-1 around the phosphorylation site of Ser259

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**In skeletal muscle, isoform 1 is more abundant than isoform 2.