
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse Phospho-SHB (Tyr246) Polyclonal Antibody

多克隆抗体

产品货号: PA01629

第三版

描述:SHB (SH2 Domain Containing Adaptor Protein B) is a Protein Coding gene. Diseases associated with SHB include Sulfhemoglobinemia. Among its related pathways are Development VEGF signaling via VEGFR2 - generic cascades and RET signaling. GO annotations related to this gene include SH3/SH2 adaptor activity. An important paralog of this gene is ENSG00000255872.

抗原:Synthesized peptide derived from human Shb around the phosphorylation site of Tyr246

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Widely expressed.