

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse Phospho-EIF4B (Ser422) Polyclonal Antibody

多克隆抗体

产品货号: PA01646

第三版

描述: EIF4B (eukaryotic translation initiation factor 4B) is a protein-coding gene. GO annotations related to this gene include RNA binding and translation initiation factor activity. An important paralog of this gene is EIF4H.

抗原: Synthesized peptide derived from human eIF4B around the phosphorylation site of Ser422

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: