

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-Ephrin B1/2 (Tyr330) Polyclonal Antibody

多克隆抗体

产品货号: PA01650

第三版

描述:EFNB1 (Ephrin B1) is a Protein Coding gene. Diseases associated with EFNB1 include Craniofrontonasal Dysplasia and Craniosynostosis. Among its related pathways are Developmental Biology and EphB-EphrinB Signaling. GO annotations related to this gene include ephrin receptor binding. An important paralog of this gene is EFNB2. EFNB2 (Ephrin B2) is a Protein Coding gene. Diseases associated with EFNB2 include Nodular Regenerative Hyperplasia and Intracranial Hemorrhage In Brain Cerebrovascular Malformations. Among its related pathways are Developmental Biology and Spinal Cord Injury. GO annotations related to this gene include receptor binding and ephrin receptor binding. An important paralog of this gene is EFNB1.

抗原:Synthesized peptide derived from human Ephrin-B1/2 around the phosphorylation site of Tyr330

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Heart, placenta, lung, liver, skeletal muscle, kidney, pancreas.Lung and kidney.