

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-FOXO1/3 (Ser322/S325) Polyclonal Antibody

多克隆抗体

产品货号: PA01660

第三版

**描述:**FOXO1 (Forkhead Box O1) is a Protein Coding gene. Diseases associated with FOXO1 include Rhabdomyosarcoma 2, Alveolar and Rhabdomyosarcoma. Among its related pathways are Constitutive Signaling by AKT1 E17K in Cancer and Common Cytokine Receptor Gamma-Chain Family Signaling Pathways. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and chromatin binding. An important paralog of this gene is FOXO3. FOXO3 (Forkhead Box O3) is a Protein Coding gene. Diseases associated with FOXO3 include Rhabdomyosarcoma and Lung Cancer. Among its related pathways are Development IGF-1 receptor signaling and Constitutive Signaling by AKT1 E17K in Cancer. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and protein kinase binding. An important paralog of this gene is FOXO1.

**抗原:**Synthesized peptide derived from human FoxO1/3 around the phosphorylation site of Ser322/S325

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**