
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-ICK (Tyr159) Polyclonal Antibody

多克隆抗体

产品货号: PA01699

第三版

描述:ICK may play a key role in the development of multiple organ systems and particularly in cardiac development.

抗原:Synthesized peptide derived from human Intestinal Cell Kinase around the phosphorylation site of Tyr159

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart, brain, placenta, pancreas, thymus, prostate, testis, ovary, small intestine and colon, with highest levels in placenta and testis. Not detected in spleen. Also expressed in many cancer cell lines.