

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-RAF1 (Ser301) Polyclonal Antibody

多克隆抗体

产品货号: PA01702

第三版

描述:Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. Part of the Ras-dependent signaling pathway from receptors to the nucleus. Protects cells from apoptosis mediated by STK3.

抗原:Synthesized peptide derived from human Raf-1 around the phosphorylation site of Ser301

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:In skeletal muscle, isoform 1 is more abundant than isoform 2.